

## Hill-Rom Affinity IV P3700 Birthing Bed Spec Sheet



### **DIMENSIONS:**

BED LENGTH – FULL.....	90” (229 cm)
BED WIDTH – SIDERAILS UP.....	39” (99 cm)
BED WIDTH – SIDERAILS DOWN.....	36” (91 cm)
HEIGHT – LOW (WITH MATTRESS).....	23” (58 cm)
HEIGHT – HIGH (WITH MATTRESS).....	39” (99 cm)
HEIGHT – MAXIMUM OF SEAT SECTION IN TREND-LIKE POSITION.....	40” (102 cm)
MATTRESS – LENGTH.....	78” (198 cm)
MATTRESS – WIDTH.....	35” (89 cm)
MATTRESS – THICKNESS (HEAD SECTION).....	5” (12.7 cm)
MATTRESS – THICKNESS (FOOT SECTION).....	3.5” (8.9 cm)
CRITICAL ANGLES – MAXIMUM HEAD ELEVATION.....	63°
CRITICAL ANGLES – MAXIMUM SEAT ELEVATION.....	15°
CRITICAL ANGLES – MAXIMUM TREND-LIKE POSITIONING.....	8°
PATIENT WEIGHT – BED.....	500 LBS (227 kg)
PATIENT WEIGHT – FOOT SECTION.....	400 LBS (181 kg)
CASTER SIZE.....	6”(15 cm)/8”(20cm)