

KCI BariKare Bariatric Bed



DIMENSIONS:

OVERALL LENGTH.....	87 ½" (222 cm)
WIDTH – SIDERAIL IN NORMAL POSITION	43 ¼" (110 cm)
WIDTH – SIDERAIL IN EXTENDED POSITION.....	54" (137.1 cm)
WIDTH – SIDERAIL IN TRANSPORT POSITION.....	40 ¼" (102 cm)
MINIMUM HEIGHT (TO "HARD PAN" SURFACE).....	21 ¼" (54 cm)
MAXIMUM HEIGHT (TO "HARD PAN" SURFACE).....	27 ½" (70 cm)
MATTRESS.....	80" (203 cm) L
.....	36" (92 cm) W
.....	5" (13 cm) H
CASTER DIAMETER.....	5" (13 cm)
SIDERAIL HEIGHT.....	21 ½" (54.6 cm)
SIDERAIL LENGTH.....	48" (121.9 cm)
BED WEIGHT.....	approx. 665 lbs (302 kg)
MAXIMUM PATIENT WEIGHT.....	850 lbs (386 kg)
ARTICULATION HEAD FOWLER.....	65°
ARTICULATION FOOT FOWLER.....	-72°
TRENDELENBURG.....	+10°
REVERSE TRENDELENBURG.....	-12.5°
CARDIAC CHAIR - FOWLER.....	65°
CARDIAC CHAIR - TRENDELENBURG.....	+10°
CARDIAC CHAIR - REVERSE TRENDELEBURG.....	-12.5°