

Stryker Epic II 2030 Critical Care Bed

DIMENSIONS:

OVERALL LENGTH.....	91" (231 cm)
OVERALL WIDTH – SIDERAILS UP.....	42 ½" (108 cm)
OVERALL WIDTH – SIDERAILS DOWN.....	40" (102 cm)
HEIGHT – LOW (TO LITTER TOP) – STANDARD.....	18" (46 cm)
HEIGHT – HIGH (TO LITTER TOP) – STANDARD.....	32 ½" (83 cm)
KNEE/GATCH ANGLE.....	0° to 30°
BACK/FOWLER ANGLE.....	0° to 90°
TRENDELENBURG.....	-12°
REVERSE TRENDELENBURG.....	+12°
PATIENT SURFACE.....	35" x 84" (89 x 213 cm)
PATIENT WEIGHT.....	500 lbs (227 kg)