

Stryker Renaissance Series 1731 PACU Stretcher



DIMENSIONS:

OVERALL LENGTH.....	84" (213.3 cm)
OVERALL WIDTH.....	34 ½" (87.6 cm)
HEIGHT – LOW.....	21 ½" (54.6 cm)
HEIGHT – HIGH.....	36" (91.4 cm)
KNEE GATCH ANGLE.....	0° to 30°
FOWLER ANGLE.....	0° to 90°
TRENDELENBURG.....	+18°
REVERSE TRENDELENBURG.....	-18°
PATIENT WEIGHT.....	500 lbs (227 kg)