

## Hill-Rom Birthing Bed

### Model: Affinity I, II

#### Feature

Overall length  
 Overall height – headboard – high position  
 Overall height – low position (without mattress)  
 Overall height - high position (without mattress)  
 Total hilow travel  
 Overall width (siderails down)  
 Siderail length  
 Clearance inside siderails  
 Siderail height (without mattress)  
 Mattress size  
 Caster wheel  
 Wheel base, center-to-center - length  
 Wheel base, center-to-center - width – head end  
 Wheel base, center-to-center – width – foot end  
 Clearance under base - high position  
 Clearance under base - low position  
 Trendelenburg angle inclination  
 Maximum patient weight  
 Foot section, maximum lift weight  
 Head section, maximum lift weight

#### Dimension

90" (228.6 cm)  
 53 ¼" (136 cm)  
 18 ½" (47 cm)  
 35 ¾" (91.44 cm)  
 17 ¼" (44.45 cm)  
 34" (86.4 cm)  
 28" (71 cm)  
 38" (96.5 cm)  
 13" (33 cm)  
 34" x 78" x 4" (86.4 cm x 198 cm x 10 cm)  
 6" (15.2 cm)  
 50" (127 cm)  
 24 ½" (64.5 cm)  
 29 ½" (75 cm)  
 5" (12.7 cm)  
 2 ½" (6.35 cm)  
 8°  
 500 lbs (227 kg)  
 400 lbs (181.6 kg)  
 200 lbs (91 kg)