

Hill-Rom TotalCare ICU Bed

Feature

Height – From Floor to Top of Articulating Deck
(bottom of mattress) in Low Position

Height – From Floor to Top of Articulating Deck
(bottom of mattress) in High Position

Minimum Under Bed Clearance

Minimum clearance between articulating deck foot section
and the floor for any bed position, except the tilt-table and the
chair egress positions with the foot section at 90°

Overall Width – Siderails Raised

Overall Width – Siderails Lowered

Overall Width – Distance between Casters

Overall Length

Overall Length – Roller Bumpers Installed

Overall Length – Transport Shelf Installed

Overall Foot Section Width

Fully Extended Length of the Foot Section
from pivot to inside of the footboard

Articulating Deck Length

Maximum Head Elevation

Maximum Knee Elevation

Maximum Foot Elevation

Maximum Foot Retract

Maximum Trendelenburg

Maximum Reverse Trendelenburg

Dimension

Between stepped sides,
14.5-15.5" (36.8-39.4 cm)

92" (233.7 cm)

1.75" (4.45 cm)

2.0" (5.1 cm)

40" (101.6 cm)

36.5" (92.7 cm)

25.5" (64.8 cm)

91.5" (232.4 cm)

92.5" (235.0 cm)

93.5" (237.5 cm)

28.88-29.13" (73.36-73.99 cm)

26.75-27.25" (67.95-69.22 cm)

81.75-82.25" (67.95-69.22 cm)

75°

20°

85° down

12" (30 cm)

12°

12°

Hill-Rom TotalCare ICU Bed

Maximum Emergency Trendelenburg	20°
Bed Safe Working Load (patient and equipment)	500 lb (226 kg)
Bed Weight with Minimum Options	525 lb (238 kg)
Bed Weight with Maximum Options	566 lb (257 kg)
Short Stay Surface Weight	36 lb (16 kg)
Treatment Surface Weight	25 lb (11 kg)
Casters	5" (13 cm) diameter
Sleep Surface Length	83" (211 cm) – primary section 56" (142 cm), and foot section 27" (69 cm)
Sleep Surface Width	Primary 35" (89 cm) and foot 31" (79 cm)
Primary Sleep Surface Thickness	9" (23 cm) at center, and 8" (20 cm) at edges (including step in the deck)
Sleep Surface Height	5" (13 cm)
Maximum Scale Weight	400 lb (181 kg)



Piedmont Medical Incorporated

500 W. Kapp St.

Dobson, NC 27017

www.piedmontmedicalinc.com

Phone: 800-433-3255

Fax: 877-433-3255