

MODEL 1060

GENERAL DIMENSIONS

Height range (to litter top)	High: 35.5"	Low: 22.0"
Caster diameter	8.0"	
Articulation		
Fowler	0° - 90°	
Trend./rev. Trend.	±12°	
Maximum weight capacity	500 lbs.	
Overall width	31.75"	
Overall length	87.5"	
Patient surface width		
Torso mattress	26.0"	
Foot-end mattress	30.0"	
Patient surface length	79.0"	
Fold-down siderail height	13.0"	
Fold-down siderail length	34.25"	
Transfer system width	8.0"	
Transfer system length	32.0"	

Stryker reserves the right to change specifications without notice.



STANDARD FEATURES

- ▶ 8" Omni Surface™ casters
- ▶ 3" T-Shaped mattress
- ▶ 26" Wide patient surface
- ▶ Steel pneumatic fowler/removable foot end litter
- ▶ Dual side-mounted foot control hydraulics
- ▶ Fifth-wheel steering
- ▶ Fluid collection basin
- ▶ Four IV receptacles
- ▶ Four spindle fold-down siderails
- ▶ Four-way adjustable foot rests
- ▶ Four-wheel ring brake system with dual end activators
- ▶ Hands-free Trendelenburg
- ▶ Integrated patient transfer boards
- ▶ Integrated storage tray with O₂ bottle holder
- ▶ Pop-up steering handles
- ▶ Roller bumpers

OPTIONAL FEATURES

- ▶ 4" Enhanced Comfort Mattress
- ▶ 3" or 4" Mattress
- ▶ 26" Patient width
- ▶ C-Spine lateral cassette film holder
- ▶ Fluid collection basin
- ▶ IV poles
 - permanently attached two- or three-stage Infusion II™ folding IV pole
 - removable IV pole
- ▶ IV caddy
- ▶ Leg supports
- ▶ Locking caster steering
- ▶ Procto board
- ▶ Patient restraints
- ▶ Serving/instrument tray
- ▶ Upright O₂ bottle holder

LABEL OPTIONS

Ambulatory Surgery	Maternity	Cath Lab
Nuclear Medicine	Emergency	PACU
Endoscopy	Radiology	Extended Stay
Surgery	GI Lab	Transport