

MODEL 5000 Series

GENERAL DIMENSIONS

Overall bed dimensions	Width: 40.0"	Length: 91.0"
Low height (top of seat litter)	17.0"	
High height (top of seat litter)	35.0"	
Mattress dimensions	Width: 33.0"	Length: 84.0"
Mattress depth	Bed: 5.0"	Foot: 3.0"
Seat depth	10.0" - 14.0"	
Caster diameter	6.0"	
Critical angles		
Head elevation	70°	
Pelvic tilt	8°	
Maximum Trendelenburg	8°	
Lift capacity		
Bed	500 lbs.	
Foot section	300 lbs.	
Head section	250 lbs.	

Stryker reserves the right to change specifications without notice.



STANDARD FEATURES

- 3" Thick foot mattress
- 5" Thick enhanced comfort mattress
- 6" Omni Surface™ casters
- 17" Low bed height
- 33" Bed width with four motor functions (head/bed/foot/Trendelenburg)
- Fluid basin
- Four-wheel ring brake system
- Glideaway foot supports with 6" adjustment and six-way adjustability
- Lightweight, self-standing foot section
- In-rail motion controls
- Locking caster steering
- Manual backup
- Night-light
- On-command pelvic tilt
- Powered leg support capability
- Removable headboard with wood laminate accent option
- Self-storing handgrips
- Single-button Trendelenburg
- Three-position siderails

OPTIONAL FEATURES

- 8" Omni Surface™ casters with wheel covers
- 500XL foot section
- Armboard
- Attachable/storable leg rests
- Basin liners
- CPR board
- CPR release
- Labor bar
- Leg supports
- Linen for head and foot mattresses
- Lumbar and perineal pillows
- Motion control pendant
- Multipurpose shelf
- Permanent folding IV pole
- Siderail communications:
Nurse call with speaker
Nurse call, TV/radio, volume, speaker, room/reading & control lock-out
- Stainless steel basin
- Translating fowler to adjust seat depth from 14" to 10"
- Wood headboard

ELECTRICAL

- 120 VAC, 60 Hz, 4 Amp: 105 to 135 VAC operating range
Optional: 220 VAX, 60 Hz, 6.5 Amp: 210 to 240 VAC operating range
- CUL
- Current leakage less than 100 microamperes
- Hospital grade plug and 3-wire heavy duty cord
- IEC 801, EMC compatible
- Four motor functions: bed, head, foot and Trendelenburg
- UL listed