

Description	Specification
Height—from floor to top of support surface (6" casters, high bed position)	39.3" (998 mm)
Height—from floor to top of support surface (6" casters, low bed position)	21.3" (541 mm)
Height—from floor to top of support surface (8" casters, high bed position)	40.0" (1016 mm)
Height—from floor to top of support surface (8" casters, low bed position)	22.5" (571 mm)
Overall length (without roller bumpers)	90" (2286 mm)
Overall length (with roller bumpers)	94" (2388 mm)

Description	Specification
Overall width (siderails in stored position)	34" (864 mm)
Overall width (siderails in up position)	38.5" (978 mm)
Hilow range of travel	17.5" (444 mm)
Sleep surface width and length	34" x 80" (864mm x 2032 mm)
Head elevation—maximum	57°
Knee elevation—maximum	32°
Trendelenburg—maximum	12°
Reverse Trendelenburg—maximum	12°
Patient weight—maximum (standard applications)	500 lbs
Safe working load—maximum (European applications)	(180 kg)
Casters	8" diameter x 1.5" wide (203 mm diameter x 38 mm wide)
Casters (optional)	6" diameter x 1.5" wide (152 mm diameter x 38 mm wide)

<b>Model</b>	<b>Description</b>
892	Bed without scale and air system
894	Bed with scale only
895	Bed with air system only
896	Bed with scale and air system

<b>Siderail Option</b>	<b>Module Number</b>	<b>Description</b>
50	380	Bed functions
51	381	Bed functions, and nurse call
52	382	Bed functions, nurse call, and entertainment
53	383	Bed functions, nurse call, entertainment, and lighting
59	389	Bed functions—upgradable to nurse call, entertainment, and lighting
60	380	Bed functions, and isolation transformer
63	383	Bed functions, nurse call, entertainment, lighting, and isolation transformer