

Reconditioned Stretchers



Stryker 1711



Stryker 721



HillRom GPS



MidMark 535

Piedmont Medical invites you to take a look at the value-added products, services and innovative solutions we have created for the healthcare industry. We offer a broad range of product options which provide our customers with excellent flexibility. At the same time, we remain committed to keeping on-hand inventory levels available for immediate purchase. These attributes, along with the company's attention to detail and quality of workmanship, have made our product easy to buy and simple to use.



AM 1010H
Overbed Table
w/vanity (New)



PMI20010
Bedside Cabinet
(New)



C1061-10
Guest Chair
(New)



500 Kapp St.
Dobson, NC 27017
(800) 433-3255
Fax: (877) 433-3255

www.piedmontmedicalinc.com
info@piedmontmedicalinc.com



Mission Statement

“Ensuring that our customers receive **VALUE** through our Commitment to **Quality Products and Efficient Core Business Principles** while effectively responding to opportunities with a sense of urgency.”

Founded in 1992, Piedmont Medical, Inc. is one of the industry’s leading suppliers of reconditioned, Hill-Rom and Stryker hospital beds, stretchers, ICU and birthing beds. Our reconditioning process, formulated to meet and exceed the demands of many medical environments, is preferred by medical professionals and consumers alike for its quality and durability.

Discover the value of our product. We invite you to purchase with confidence and pride.

ROOM SERVICE



It's as easy as “one two three”. Select product from our Room Service menu with our standard finish options and we will give you a “one-room” price. The six piece group, as pictured above, includes the following: reconditioned Hill-Rom 850 bed, new #500 mattress, AM4828U overbed table, PMI20010 bedside cabinet, C540-40 high-back chair and a C51-5 guest chair.

Redefining the Meaning of “Reconditioned”

All beds are visually inspected for bent or broken parts. Any bent or broken parts are repaired or replaced so that operating features function properly.

Our beds are completely disassembled.

Our beds are cleaned inside and outside of the frame and all stickers or decals are removed.

Our bed frames, arm rails and pans (if applicable) are sanded completely.

All beds and parts are then cleaned again to remove any dust or dirt.

Bed frames, arm rails and pans are painted in our professional paint booth. All parts and frames dry for 24 hours.

Beds are then reassembled and include the following NEW parts:

- NEW Tente casters
- NEW plastic parts on arm rails
- NEW head and foot boards
- NEW labels or decals
- NEW capacitor and lubricant for motors

A final inspection of the beds operating features is performed and it is tested for electrical safety to ensure performance meets manufacturers’ standards.

Beds are blanket wrapped, covered in protective plastic and positioned in a custom stand for stability during shipment. Head and foot boards are shipped unattached and seperately.



LOOKS LIKE NEW, WORKS LIKE NEW OR IT NEVER GETS TO YOU!

