

Stryker 1710/1711 Stretcher

<u>Feature</u>	<u>Dimension</u>
Maximum Weight Capacity	500 pounds
Overall Bed Length	84" (213.4 cm)
Overall Bed Width	34.5" (87.6 cm)
Minimum Bed Height	22" (55.9 cm)
Maximum Bed Height	35.5" (90.2 cm)
Fowler Angle	0 to 90°
Knee Gatch Angle	0 to 30°
Trendelenburg / Reverse Trendelenburg	+18° to -18°