

Beds, Sleep Surfaces, Stretchers and Furniture  
800-STRYKER (787-9537) or 269-329-2100  
www.strykermedical.com



MODEL 721

## FEATURES AND BENEFITS

- Optional, built-in transfer board (long or short available)
- Fifth wheel activates for easier traction and cornering
- Permanently-mounted IV caddy holds freestanding IV poles
- Four-wheel, steel-ring brakes with optional four-sided controls

### COLOR OPTIONS



MODEL 720 (FIXED HEIGHT)

# TRANSPORT STRETCHERS

## STANDARD FEATURES

- 3" (8 cm) mattress
- 8" (20 cm) Omni Surface™ casters
- 24" (61 cm) wide patient surface
- Fold-to-foot siderails
- Four IV receptacles
- Four-wheel, steel-ring brakes with dual end activators
- Integrated storage tray with oxygen bottle holder
- Roller bumpers
- Stationary foot-end litter

## MODEL 721: ALL STANDARD FEATURES PLUS

- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Dual side-mounted foot control hydraulics
- Hands-free Trendelenburg

## MODEL 720: ALL STANDARD FEATURES PLUS

- Dual pedestal design
- Fixed height of 28", 30" or 32" (71, 76 or 81 cm)
- Manual Trendelenburg

## OPTIONAL FEATURES

- BIG WHEEL® (721 only)
- 3" (8 cm) or 4" (10 cm) Enhanced Comfort Mattress
- 4" (10 cm) mattress
- 26" (66 cm) wide patient surface
- C-spine cassette film holder
- Constant descent hydraulics (721 only)
- Defibrillator tray
- Defibrillator tray/footboard/extender
- Footboard
- Footboard/chartholder
- Four-sided brake & steer controls
- Integrated patient transfer boards (long or short)
- IV caddy
- IV poles
  - permanent folding two- or three-stage Infusion™ II
  - tethered
  - standard removable
- Patient restraints
- Pneumatic backrest (720 only)
- Retractable fifth-wheel steering
- Serving/instrument tray
- Siderail pads
- Upright oxygen bottle holder
- Wheel covers
- Extended warranty available


## SPECIFICATIONS

• <b>Model Number</b>	721		720	
• <b>Overall Length</b>	83"	(211 cm)	81"	(206 cm)
• <b>Overall Width</b>	30"	(76 cm)	30"	(76 cm)
• <b>Weight Capacity</b>	500 lbs	(228 kg)	500 lbs	(228 kg)
• <b>Height Range (to litter top)</b>				
<i>High</i>	36"	(91 cm)	N/A	
<i>LOW (with short transfer board)</i>	21.5"	(55 cm)	N/A	
<i>LOW (with long transfer board)</i>	24.5"	(62 cm)	N/A	
• <b>Height</b>				
<i>Fixed</i>	N/A		28", 30" or 32"	(71, 76, 81 cm)
• <b>Litter Positioning</b>				
<i>Backrest</i>	0° - 90°		0° - 90°	
<i>Trend./Reverse Trend.</i>	±18°		12° Trend.	
• <b>Patient Surface</b>				
<i>Standard</i>	24" x 76"	(61 x 193 cm)	24" x 76"	(61 x 193 cm)
<i>Optional</i>	26" x 76"	(66 x 193 cm)	26" x 76"	(66 x 193 cm)
• <b>Siderails</b>	13" x 55"	(33 x 140 cm)	13" x 55"	(33 x 140 cm)
• <b>Transfer Board</b>				
<i>Long</i>	10" x 64.5"	(25 x 164 cm)	10" x 64.5"	(25 x 164 cm)
<i>Short</i>	8" x 44"	(20 x 112 cm)	8" x 44"	(20 x 112 cm)

## LABEL OPTIONS

Ambulatory Surgery	Maternity	Cath Lab	Nuclear Medicine
Emergency	PACU	Extended Stay	Radiology
Surgery	Trauma	Cardiology	Transport
Same Day Surgery	UltraSound	GI Lab	Outpatient Surgery
Recovery	Endoscopy		

*Stryker reserves the right to change specifications without notice.*



**stryker**®

For Beds, Sleep Surfaces, Stretchers and Furniture

6300 South Sprinkle Road, Kalamazoo, Mich. 49001 USA  
 269.329.2100 800.STRYKER 800.FAX.STRY  
 www.strykermedical.com

© Stryker Corporation SCG-15M-0903 LSS-355-10B