

## Hill-Rom TranStar® Stretchers



### At a glance

- Low height improves on/off transfer and increases patient safety
- One-Step® TuckAway siderails for zero-gap transfers (Standard on PC-300 through PC-550 value packages)
- Steering Plus™ 5th wheel system for superior maneuverability and sharp cornering (PC-100/PC-200 value packages feature 4-corner brake and steer controls)
- Ergonomic push handles to reduce caregiver strain during transport (Standard on PC-200 through PC-550 value packages)
- GentleRide™ shock absorbing frame and siderail system for improved patient comfort
- TriPedal foot controls allow even lowering of patients (Standard on PC-300 through PC-550 value packages)
- Permanent IV pole (Standard on PC-200 through PC-550 value packages)

Designed with patients and caregivers in mind, TranStar stretchers offer such features as an advanced steering system, ergonomic push handles, a shock absorbing frame and siderail system and an array of options and accessories. Coupled with the support and service you can only get from Hill-Rom, TranStar stretchers are the ideal solution for patient care and transport throughout the healthcare facility.



### Procedural stretcher: Effortless mobility and high durability (PC-300 through PC-550 value packages)

- Ideal for use in emergency, OR, PACU and outpatient.
- Siderails quickly tuck away for zero-gap transfers.
- Ergonomically correct transport handles reduce caregiver strain during patient transport.
- Steering Plus™ 5th wheel system for true tracking and easy, sharp cornering.
- Low height (21.5") improves on/off transfers and increases patient safety.
- Tri-pedal foot controls allow even lowering of patients.

### Procedural value stretcher: The economical choice (PC-100/PC-200 value packages)

- Four-corner brake/steer controls offer improved maneuverability.
- Ergonomic transport handles reduce caregiver strain during transport (PC-200 value package only).
- Siderails designed to minimize pinched fingers and tubes.
- Easily placed into emergency Trendelenburg position by lifting the foot end of the stretcher.

# Hill-Rom TranStar® Stretchers



Optional fiber optic exam light



## OB/GYN stretcher: Easy patient positioning

- Ideal for use in emergency and maternity.
- Store-away foot section eliminates heavy lifting and storage concerns.
- Sliding litter enables quick and easy patient positioning.
- Integrated foot supports and catch basin offer convenience and time savings.
- Optional fiber optic exam light easily positionable for precision and convenience.




Full length standard X-ray cassette platform under sleep deck



## Trauma stretcher: Maximum patient access

- Ideal for use in trauma and emergency.
- X-ray capabilities from head-to-toe without moving the patient, saving time and reducing injury risk.
- Full-length platform accepts all standard X-ray cassettes and grids.
- X-ray cassette holder for upright chest positions allows for increased efficiency.
- Drainage bags/chest tube suction mounts at all four corners.



Siderail controls to raise and lower the head and knee sections are easy to read and activate



Choice of superior or temporal wrist rest with single-piece flat or two-piece flat/concave cushion (superior wrist rest with two-piece flat/concave cushion shown)

**Electric stretcher: Comfort and convenience of a hospital bed, mobility of a stretcher**

- Ideal for use in PACU, ER observation, and outpatient.
- Siderail controls allow patients to adjust the stretcher for their individual comfort.
- Staff control panel at footend allows for greater caregiver control.

**Surgical stretcher: Precise surgeon functionality**

- Ideal for use in outpatient surgery, OR, and freestanding surgery centers.
- Articulating head section and wrist rest options to suit a variety of surgeon preferences.
- Surgical brake/steer pedal and twin wheel casters increase stability and provide maximum foot space.
- Surgical rail allows attachment of gas delivery/drape support devices, PACU extender/armboards and surgical hand tables.

# Hill-Rom TranStar® Stretchers



## TranStar® Stretcher value package configurations

Model	26" width	30" width	Ergonomic Push Handles	IV Pole(s)	5th Wheel Steering	Knee Gatch	Radio-lucent Surface	Footend Pedals	Factory Mounted Light	3" Thick	Mattress Options			30" Width
											4" Thick	5" Thick	26" Width	
PC-550 Procedural		x	x	2	x	x	x	x		Std	Op	Op		Std
PC-500 Procedural	x		x	2	x	x	x	x		Std	Op	Op		Std
PC-450 Procedural		x	x	1	x	x				Std	Op	Op		Std
PC-400 Procedural	x		x	1	x	x				Std	Op	Op		Std
PC-350 Procedural		x	x	1	x					Std	Op	Op		Std
PC-300 Procedural	x		x	1	x					Std	Op	Op		Std
PC-200 Procedural	x		x	1						Std	Op	Op		Std
PC-100 Procedural	x									Std	Op	Op		Std
TR-500 Trauma	x		x	1	x		x	x		Std	Op	Op		Std
OB-550 OB/GYN		x	x	1	x		x		x	Std				Std
OB-450 OB/GYN		x	x	1	x		x			Std				Std
SG-500 Surgical	x		x	1	x	x					Std			Std
EL-550 Electric		x	x	1	x	x					Std	Op		Std



### Accessories

Telescopic IV Pole  
ISS Transfer Pole  
Socket Adapter (for ISS Transfer Pole)  
IV Transporter  
Upright Chest X-ray Cassette Holder (Trauma and OB/GYN)  
Lateral X-ray Cassette Holder  
Liquid Oxygen Tank Holder  
Vertical Oxygen Tank Holder  
Chart Holder  
Paper Roll Dispenser (standard/wide width)

Transport Straps (1 strap per package)  
3" and 4" Foot Extenders  
Siderail Pads  
Mattress Pad  
Stretcher Pillow  
Patient Tray  
Utility Tray (standard/wide width)  
Detachable Head/Footboard (headboard for OB/GYN)  
Convertible Footboard  
Armboard

Ankle Stirrups  
Exam Light Upgrade (OB/GYN only)  
Detachable Catch/Fluid Basin (OB/GYN only)  
Telescoping Calf Supports (OB/GYN only)  
Temporal Wrist Rest (Surgical Only)  
Superior Wrist Rest (Surgical Only)  
Gas Delivery/Drape Support Device (Surgical Only)  
Head Restrain Strap (Surgical Only)  
Integrated PACU Extenders (Surgical Only)

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2003 Hill-Rom Services, Inc.  
ALL RIGHTS RESERVED.

CTG470rb 6/13/03

Note: CAL 129 Flammability requirement is an option on all mattresses.

Labeling: Standard decals include: Emergency, Endoscopy, Cath Lab, GI Lab, PACU, Transport, Maternity, Radiology, Surgery, and Ambulatory Surgery

Custom decals: See your Hill-Rom Customer Service Agent or Account/Territory Manager if custom decals are desired.