

<b>Feature</b>	<b>Dimension</b>
Overall length	90" (228.6 cm)
Overall height - headboard - high position	53 1/4" (136 cm)
Overall height - low position (without mattress)	18 1/2" (47 cm)
Overall height - high position (without mattress)	35 3/4" (91.44 cm)
Total hilow travel	17 1/4" (44.45 cm)
Overall width (siderails down)	34" (86.4 cm)
Siderail length	28" (71 cm)
Clearance inside siderails	38" (96.5 cm)
Siderail height (without mattress)	13" (33 cm)
Mattress size	34" x 78" x 4" (86.4 cm x 198 cm x 10 cm)
Caster wheel	6" (15.2 cm)
Wheel base, center-to-center - length	50" (127 cm)
Wheel base, center-to-center - width - head end	24.5" (64.5 cm)
Wheel base, center-to-center - width - foot end	29.5" (75 cm)
Clearance under base - high position	5" (12.7 cm)
Clearance under base - low position	2.5" (6.35 cm)
Trendelenburg angle inclination	8°
Maximum patient weight	500 lb (227 kg)
Foot section, maximum lift weight	400 lb (181.6 kg)
Head section, maximum lift weight	200 lb (91 kg)