


Standard Comfort Stretcher Pads

Standard Comfort

Standard Comfort Stretcher Pads are available in seven 'in stock' sizes with universal corners - 2radius/2angle. Medium firm high density comfort foam covered with electrically conductive vinyl provides a comfortable latex free, bacteria and fluid resistant pad for all gurney and transport stretcher applications.

Upgrade Option: request 4-way stretch vinyl fabric to provide extra comfort for your patients.

Custom: Custom stretcher pads are made to your specifications, choice of foam, corner configurations, vinyl covering and thicknesses. Available in both Standard and Kova Pressure Management comfort levels.




MOST POPULAR SIZE PADS:		
MAKE	MODEL	Dimensions
Standard Stretcher Pads	SP-5000	24" x 76" x 3"
	SP-5004	24" x 75" x 4"
	SP-5017	24" x 77" x 4"
	SP-5051	27" x 75" x 4"
	SP-5054	27" x 77" x 4"
	SP-5057	29" x 75" x 4"
	SP-5060	30" x 76" x 4"

*Available for all makes and models of stretchers.
Most popular styles shown in table above, are in stock.*



Angle Corners



Radius Corners