

Long Term Care

**3 Layered Mattress with
Visco Foam Heel Section
and a 7 Year Warranty**



COMPARE TO

*Direct Supply® Panacea® Tender
Invacare® Soft Forms
Span America® Gromattress
Blue Chip® Pre-vent II Series

accessories

15770C	Replacement Cover	1/ea
15777C	Replacement Cover	1/ea
15784C	Replacement Cover	1/ea
15785C	Replacement Cover	1/ea
15876C	Replacement Cover	1/ea
15877C	Replacement Cover	1/ea



Gravity 7 Long Term Care Pressure Redistribution Mattress

Item #	Raised Side Rails	Product Description	Product Dimensions	Product HCPCS	UOM
15876	No	Multi-Layered/Multi-zoned Foam Mattress	36" (W) x 76" (L) x 6" (H)	E0184	1/ea
15770	No	Multi-Layered/Multi-zoned Foam Mattress	36" (W) x 80" (L) x 6" (H)	E0184	1/ea
15784	No	Multi-Layered/Multi-zoned Foam Mattress	36" (W) x 84" (L) x 6" (H)	E0184	1/ea
15877	Yes	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out	36" (W) x 76" (L) x 6" (H)	E0184	1/ea
15777	Yes	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out	36" (W) x 80" (L) x 6" (H)	E0184	1/ea
15785	Yes	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out	36" (W) x 84" (L) x 6" (H)	E0184	1/ea

*Direct Supply® and Panacea® are registered trademarks of Direct Supply, Inc.

specifications

Cover Material	Stretch, Fluid-Resistant Nylon Top Waterproof Vinyl Bottom
Mattress Layers of Foam	3
Side Rails	Yes
Dedicated Heel Section	Yes
Fire Rating	CFR 16 Part 1633 Federal Fire Code
Weight Capacity	350 lbs.
Warranty	7 Year Limited Foam 1 Year Cover

- Drive's most economical yet efficient therapeutic foam mattress maximizes value and comfort with 3 full layers of foam
- Top layer uses unique die-cutting and a true memory foam heel section to optimize zoned pressure redistribution
- Bottom layer provides full length horizontally scored articulation cuts extending the durability and life of the mattress significantly
- Stretch cover is fire retardant, fluid and tear resistant, low shear and easy to clean
- High Density Foam Perimeter provides solid edge support and facilitates safe patient transfers
- Available with or without raised side rails