


# Long Term Care

**15 Year Warranty**

## COMPARE TO

\*Direct Supply® Panacea® Clinical  
Invacare® Therapy 3000  
Span America® GeoMatt UltraMax  
Joerns® Prevamatt Elite


**Full Length Open-Cell  
Visco Top Layer Heel Slope  
for Top of the Line  
Prevention & Therapy**



## accessories

15970C	Replacement Cover for 15970	1/ea
15977C	Replacement Cover for 15977	1/ea
15984C	Replacement Cover for 15984	1/ea
15985C	Replacement Cover for 15985	1/ea
15996C	Replacement Cover for 15996	1/ea
15997C	Replacement Cover for 15997	1/ea

## Heel Slope



## Gravity 9 Premium Long Term Care Pressure Redistribution Mattress

Item #	Raised Side Rails	Product Description	Product Dimensions	Product HCPCS	UOM
15996	No	Multi-Layered/Multi-zoned Foam Mattress	36" (W) X 76" (L) X 6" (H)	E0184	1/ea
15970	No	Multi-Layered/Multi-zoned Foam Mattress	36" (W) X 80" (L) X 6" (H)	E0184	1/ea
15984	No	Multi-Layered/Multi-zoned Foam Mattress	36" (W) X 84" (L) X 6" (H)	E0184	1/ea
15997	Yes	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out	36" (W) X 76" (L) X 6" (H)	E0184	1/ea
15977	Yes	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out	36" (W) X 80" (L) X 6" (H)	E0184	1/ea
15985	Yes	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out	36" (W) X 84" (L) X 6" (H)	E0184	1/ea

- The ultimate in the industry for multi-layered die-cut foam and true memory foam, available with or without raised side rails
- Mattress includes a tapered heel slope and channel cut center layer of foam under the entire lower leg to offload heels and prevent heel pressure ulcers
- Top layer is full-length open-cell visco-elastic memory foam
- Bottom layer provides full length horizontally scored articulation cuts extending the durability and life of the mattress significantly
- Stretch cover is fire retardant, fluid and tear resistant, low shear and easy to clean
- High Density Foam Perimeter provides solid edge support and facilitates safe patient transfers

## specifications

Cover Material	Stretch, Fluid-Resistant Nylon Top   Waterproof Vinyl Bottom
Mattress Layers of Foam	4
Side Rails	Yes
Dedicated Heel Section	Yes

Fire Rating	CFR 16 Part 1633 Federal Fire Code
Weight Capacity	450 lbs.
Warranty	15 Year Limited Foam   1 Year Cover