

HILL-ROM 850 HOSPITAL BED ROOM PACKAGE



HOSPITAL BED (Hill-Rom 850)

FEATURES:

- Height / Floor to Top of Surface / High Position: 31" (78.7 cm)
- Height / Floor to Top of Surface / Low Position: 17" (43.2 cm)
- Overall Length / With Roller Bumpers / Fully Extended: 93-1/2" (237.5 cm)
- Overall Length / With Roller Bumpers / Fully Retracted: 79-1/2" (201.9 cm)
- Overall Length / Without Bumpers / Fully Extended: 91" (231.1 cm)
- Overall Length / Without Bumpers / Fully Retracted: 77" (195.6 cm)
- Overall Width / Siderails Up: 41-1/2" (105.4 cm)
- Total Hilow Travel: 14" (35.6 cm)
- Sleeping Surface: 36" x 82" (91.4 cm) x 208.3 cm)
- Weight / Bed with Casters: 346 lb (156.9 kg)
- Maximum Patient Weight: 500 lb (226.8 kg)

MATTRESS (Gravity 7 Long Term Care Pressure Redistribution Mattress)

- Most economical yet efficient therapeutic foam mattress maximizes value and comfort with 3 full layers of foam
- Top layer uses unique die-cutting and a true memory foam heel section to optimize zoned pressure redistribution.
- Bottom layer provides full length horizontally scored articulation cuts extending the durability and life of the mattress significantly.
- Stretch cover is fire retardant, fluid, and tear resistant, low shear and easy to clean.
- High Density Foam Perimeter provides solid edge support and facilitates safe patient transfers.
- Available with or without raised side rails.
- Cover Material: Stretch, Fluid-Resistant Nylon Top | Waterproof Vinyl Bottom
- Fire Rating: CFR 16 Part 1633 Federal Fire Code
- Weight Capacity: 350 lbs.
- Warranty: 7 Year Limited Foam | 1 Year Cover

BEDSIDE CABINET (One Drawer/One Door)

- 24.25"W x 16.75"D x 30.75"H
- One Door/One Drawer
- Casters
- Weight: 65 lbs

OVERBED TABLE (Model 4528U)

FEATURES:

- 18"W x 32"L High Pressure laminate top
- chrome plated 11 and 12 gage welded tubular steel column and base
- 2" twin wheel casters
- 30"H - 45"H adjustable in 15 one inch increments (Optional Model: 1010H)

STANDARD FEATURES:

- Chrome base finish
- Standard laminate options

Piedmont Medical, Inc.
PO Box 1278, 500 Kapp St. Dobson, NC 27017
Phone: 800-433-3255 Fax: 877-433-3255