

HILL-ROM ADVANTA HOSPITAL BED ROOM PACKAGE



HOSPITAL BED (Hill-Rom Advanta)

FEATURES:

- Overall Length: 91" (231 cm)
- Retracted Length: 78-1/2" (199.4 cm)
- Mattress Support Width: 35" (89 cm)
- Mattress Support Length: 83" (211 cm)
- Maximum Width / Siderails Up: 42-3/4" (108.59 cm)
- Maximum Width / Siderails Stored: 36" (91 cm)
- Maximum Headboard / Height: 48" (122 cm)
- Maximum Footboard / Height: 45" (114 cm)
- Maximum Siderail / Height: 45" (114 cm)
- Bed Clearance: 4-3/4" (12.07 cm)
- Head Section Inclination / Maximum: 65°
- Knee Section Inclination / Maximum: 35°
- Automatic Contour Up / Maximum: 15°
- Automatic Contour Down: 28°
- Automatic Contour: 10°
- Hilow Range: 13-3/4" (33.66 cm)
- Trendelenburg: 11°
- Reverse Trendelenburg: 11°
- Maximum Patient Weight: 500 lb (226.8 kg)

MATTRESS (Gravity 7 Long Term Care Pressure Redistribution Mattress)

- Most economical yet efficient therapeutic foam mattress maximizes value and comfort with 3 full layers of foam
- Top layer uses unique die-cutting and a true memory foam heel section to optimize zoned pressure redistribution.
- Bottom layer provides full length horizontally scored articulation cuts extending the durability and life of the mattress significantly.
- Stretch cover is fire retardant, fluid, and tear resistant, low shear and easy to clean.
- High Density Foam Perimeter provides solid edge support and facilitates safe patient transfers.
- Available with or without raised side rails.
- Cover Material: Stretch, Fluid-Resistant Nylon Top | Waterproof Vinyl Bottom
- Fire Rating: CFR 16 Part 1633 Federal Fire Code
- Weigh Capacity: 350 lbs.
- Warranty: 7 Year Limited Foam | 1 Year Cover

BEDSIDE CABINET (One Drawer/One Door)

- 24.25"W x 16.75"D x 30.75"H
- One Door/One Drawer
- Casters
- Weight: 65 lbs

OVERBED TABLE (Model 4528U)

FEATURES:

- 18"W x 32"L High Pressure laminate top
- chrome plated 11 and 12 gage welded

STANDARD FEATURES:

- Chrome base finish
- Standard laminate options

Piedmont Medical, Inc.
PO Box 1278, 500 Kapp St. Dobson, NC 27017
Phone: 800-433-3255 Fax: 877-433-3255

HILL-ROM ADVANTA HOSPITAL BED ROOM PACKAGE

- tubular steel column and base
 - 2" twin wheel casters
 - 30"H - 45"H adjustable in 15 one inch increments
- (Optional Model: 1010H)

Piedmont Medical, Inc.
PO Box 1278, 500 Kapp St. Dobson, NC 27017
Phone: 800-433-3255 Fax: 877-433-3255