

Hill-Rom Transtar 8000 Procedural Stretcher

Product Description

- Height Overall / High Position: **35" (88.9 cm)**
- Height Overall / Low Position: **21½" (54.6 cm)**
- Overall Length: **83" (210.8 cm)**
- Overall Width / Siderails up: **32" (81.3 cm)**
- Upper Frame Width / Siderails Stored: **26½" (67.3 cm)**
- Overall Siderail Length: **47" (119.4 cm)**
- Siderail Height Above Sleep Surface: **14" (35.6 cm)**
- Mattress Size: **26" x 75" (66 cm x 190.5 cm)**
- Caster Size: **8" (20.3 cm)**
- Floor-to-Base Clearance: **3½" (8.9 cm)**
- Floor-to-Hydraulic Cylinder Clearance: **2" (5 cm)**
- Base Shroud Width: **28" (71.1 cm)**
- Maximum Back Incline: **90°**
- Maximum Thigh Section Incline: **25°**
- Weight / Stretcher with Casters: **225 – 260 lb. (102 – 118 kg)**
- Maximum Patient Weight: **500 lb. (226.8 kg)**
- Standard Features: **Back Section Inclination / Thigh Section Inclination / Foot Section Inclination / IV Pole Accommodation**

Optional Features: **NEW STRETCHER PAD SOLD SEPARATELY**

