

# Hill-Rom 8050 TranStar OBGYN Stretcher

## Product Description

- Height Overall / High Position: **35" (88.9 cm)**
- Height Overall / Low Position: **21.5" (54.6 cm)**
- Overall Length: **83" (210.8 cm)**
- Overall Width / Siderails Up: **32" (81.3 cm)**
- Upper Frame Width / Siderails Stored: **26.5" (67.3 cm)**
- Overall Siderail Length: **47" (119.4 cm)**
- Siderail Height Above Sleep Surface: **14" (35.6 cm)**
- Mattress Size: **26" x 75" (66 cm x 190.5 cm)**
- Caster Size: **8" (20.3 cm)**
- Floor-to-Base Clearance: **3½" (8.9 cm)**
- Floor-to-Hydraulic Cylinder Clearance: **2" (5 cm)**
- Base Shroud Width: **28" (71.1 cm)**
- Maximum Back Incline: **90°**
- Maximum Foot Section Incline: **90°**
- Weight / Stretcher with Casters: **225 – 290 lb. (115.7 – 131.5 kg)**
- Maximum Patient Weight: **500 lb. (226.8 kg)**
- Optional Features: **NEW STRETCHER PAD SOLD SEPARATELY**

