

## Hill-Rom P1170 CareAssist Bed

### Product Description

- Height – High Position: **32 ½" (82.6 cm)**
- Height – Low Position: **16 ¾" (42.5 cm)**
- Total Length: **100" (254 cm)**
- Maximum Width (siderails up): **40" (102 cm)**
- Maximum Width (siderails stored): **40" (102 cm)**
- Caster Size: **6" (15 cm)**
- Sleep Surface: **36" x 80" (91.4 cm x 203.2 cm)**
- Trendelenburg: **16°**
- Head Section Inclination: **65°**
- Knee Section Inclination: **20°**
- Foot Section Inclination: **-23°**
- Bed Weight with Casters: **420 lb. (191 kg)**
- Maximum Safe Working Load: **400 lb. (181 kg)**

