

Hill-Rom P3200 VersaCare Bed

Product Description

- Overall Length – Foot Section Extended: **94 ½" (240 cm)**
- Overall Length – Foot Section Retracted: **82 ½" (210 cm)**
- Maximum Width (siderails stored): **37" (94 cm)**
- Maximum Width (siderails up): **40" (102 cm)**
- Maximum Headboard Height: **39 ½" (100 cm)**
- Mattress to Siderail Height: **9 ½" (24 cm)**
- Minimum Under-Bed Clearance: **3" (8 cm) 1 ¼" (3.18 cm)**
- **IntelliDrive® Transport System beds** Wheel Base: **26 ¾" x 62" (67.95 cm x 157.48 cm)**
- Surface Dimensions: **35 ½" x 86 2/5" (90 cm x 219.5 cm)**
- Caster Size: **5" (13 cm) standard, or 6" (15 cm) optional (standard on IntelliDrive® Transport System beds)**
- Bed Lift Capacity (safe working load) (includes patient weight, mattress, IV poles, and such): **550 lb. (249 kg)**
- Total Weight without Surface: **495 lb. (224 kg)**
- Head Section Inclination (maximum): **65°**
- Knee Section Inclination (maximum): **16°**
- Foot Section Inclination (maximum): **-27°**

