

Stryker 1007 M-Series Stretcher

Product Description

- Overall Length: **84" (213.4 cm)**
- Overall Width (siderails up): **33 ½" (85.1 cm) [26" litter] or 37" (94 cm) [30" litter]**
- Overall Width (siderails down): **30 ¼" (76.8 cm) [26" litter] or 30 ¾" (78.1 cm) [30" litter]**
- Height Overall – High Position: **34" (86.4 cm)**
- Height Overall – Low Position: **20 ¾" (52.7 cm)**
- Caster Size: **8" (20.3 cm)**
- Knee Gatch: **0° – 40°**
- Fowler Angle: **0° – 90°**
- Trendelenburg/Reverse Trendelenburg: **±16°**
- Minimum Under Stretcher Clearance: **6" (15 cm) nominal 1 ¾" (4.5 cm) under the hydraulic cylinders and 5th wheel**
- Maximum Patient Weight: **700 lb. (317.5 kg)**

