

Stryker 1010 Renaissance Stretcher

Product Description

- Overall Length: **83" (210.8 cm)**
- Overall Width: **31.5" (80 cm)**
- Height Overall / High Position: **35.5" (90.2 cm)**
- Height Overall / Low Position: **22" (55.9 cm)**
- Patient Surface: **26" x 75.5" (66 cm x 191.8 cm)**
- Siderails: **15" x 55" (38.1 cm x 139.7 cm)**
- Caster Size: **8" (20.3 cm)**
- Knee Gatch: **0° – 30°**
- Backrest: **0° – 90°**
- Trendelenburg: **+12°**
- Reverse Trendelenburg: **-18°**
- Maximum Patient Weight: **500 lb. (226.8 kg)**
- Bed Weight: **280 lb. (127 kg)**
- Optional Features: **NEW STRETCHER PAD SOLD SEPARATELY**

