

Stryker 1015 Big Wheel M-Series Stretcher

Product Description

- Overall Length: **84" (213.4 cm)**
- Overall Width (siderails up): **33 ½" (85.1 cm) or 37" (94 cm) [30" litter]**
- Overall Width (siderails down): **30 ¼" (76.8 cm) or 30 ¾" (78.1 cm) [30" litter]**
- Height Overall / High Position: **34" (86.4 cm)**
- Height Overall / Low Position: **20 ¾" (52.7 cm)**
- Patient Surface: **26" x 74 ¼" (66 cm x 188.6 cm) or 30" x 74 ¼" (76.2 cm x 188.6 cm)**
- Siderails: **14" x 57 ¾" (35.6 cm x 146.7 cm)**
- Caster Size: **8" (20.3 cm)**
- Knee Gatch: **0° - 40°**
- Backrest: **0° - 90°**
- Maximum Angle: **87.5°**
- Trendelenburg/Reverse Trendelenburg: **±16°**
- Maximum Patient Weight: **700 lb. (317.5 kg)**

