

Stryker 1020 Trauma Stretcher

Product Description

- Overall Length: **83" (210.8 cm)**
- Overall Width (siderails up): **31 ½" (80 cm)**
- Overall Width (siderails down): **30" (76.2 cm)**
- Height Overall – High Position: **36" (91.4 cm)**
- Height Overall – Low Position: **21.5" (55.9 cm)**
- Patient Surface: **26" x 75" (66 cm x 190.5 cm)** or **29" x 75" (74 cm x 190.5 cm)**
- Siderails: **13" x 55" (33 cm x 139.7 cm)**
- Caster Size: **8" (20.3 cm)**
- Knee Gatch: **0° – 30°**
- Backrest: **0° – 90°**
- Trendelenburg/Reverse Trendelenburg: **±18°**
- Maximum Patient Weight: **500 lb. (226.8 kg)**

