

Stryker 1061 Gynnie OBGYN Stretcher

Product Description

- Overall Length: **82" (208.3 cm)**
- Overall Width: **31 ½" (80 cm)**
- Height Overall – High Position: **35 ½" (90.2 cm)**
- Height Overall – Low Position: **22" (55.9 cm)**
- Caster Size: **8" (20.3 cm)**
- Siderails: **13" x 39" (33 cm x 99.1 cm)**
- Patient Surface: **26" x 77" (66 cm x 195.6 cm)**
- Backrest: **0° – 90°**
- Trendelenburg/Reverse Trendelenburg: **±19°**
- Maximum Patient Weight: **500 lb. (226.8 kg)**

