

Stryker 1711 Renaissance Stretcher

Product Description

- Overall Length: **84" (213.4 cm)**
- Overall Width (siderails up): **34 ½" (87.6 cm)**
- Overall Width (siderails down): **30 ½" (77.5 cm)**
- Height Overall – High Position: **35 ½" (90.2 cm)**
- Height Overall – Low Position: **22" (55.9 cm)**
- Patient Surface: **29" x 76" (73.7 cm x 193 cm)**
- Siderails: **14" x 54" (35.6 cm x 137.2 cm)**
- Caster Size: **8" (20.3 cm)**
- Knee Gatch: **0° – 30°**
- Backrest: **0° – 90°**
- Trendelenburg/Reverse Trendelenburg: **±18°**
- Maximum Patient Weight: **500 lb. (226.8 kg)**

