

# Stryker 2040 Zoom ICU Bed

## Product Description

Overall Length: **93" (236 cm)**

- Overall Width / Siderails Up: **42.5: (108 cm)**
- Overall Width / Siderails Down: **40" (102 cm)**
- Weight Capacity: **500 lbs (227 kg)**
- Standard Height Range / To Litter Top / High: **32.5" (83 cm)**
- Standard Height Range / To Litter Top / Low: **18" (46 cm)**
- Litter Positioning / Backrest: **0 to 90 degrees**
- Litter Positioning / Knee Gatch: **0 to 30 degrees**
- Trendelenburg / Reverse Trendelenburg: **+10/-12 degrees**
- Patient Surface: **35" x 84" (89 x 213 cm)**
- Caster Diameter: **Standard 6" (15 cm)**

