

Stryker 5050 Stretcher Chair

Product Description

- Overall Length: **76" (193 cm)**
- Overall Width: **30" (76.2 cm)**
- Overall Width (litter): **27 ½" (69.8 cm)**
- Height Overall / High Position: **33" (83.8cm)**
- Height Overall / Low Position: **22" (55.9 cm)**
- Patient Surface: **24" x 74" (61 cm x 188 cm)**
- Siderails: **10" x 31" (25.4 cm x 78.7 cm)**
- Caster Size: **6" (15.2 cm)**
- Leg: **0° – 80°**
- Backrest: **0° – 90°**
- Trendelenburg/Reverse Trendelenburg: **±18°**
- Maximum Patient Weight: **400 lb. (181.4 kg)**

