

# Stryker S3 Bed

## Product Description

- Overall Length: **93" (236.2 cm)**
- Overall Width – Siderails Up: **41.5" (105.4 cm)**
- Overall Width – Siderails Down: **39.5" (100.3 cm)**
- Height Overall – High Position: **30" (76.2 cm)**
- Height Overall – Low Position – Standard with 6" Casters: **16" (40.6 cm)**
- Patient Surface: **35" x 84" (89 x 213 cm)**
- Caster Size: **6" (15 cm)**
- Maximum Patient Weight: **500 lb. (227 kg)**

## Standard Features:

- Centrally-located brake pedal
- 4-wheel, steer brake with swivel lock system
- One-touch 30-degree head of bed angle and lock
- CPR Release and removable headboard
- Trend. / Reverse Trend. Controls
- Night Light
- Dual pedestal lift design
- Four IV pole and traction equipment sockets
- Three-sided bed motion controls

