

Stryker Secure II 3002

The Stryker 3002 Secure II Bed is safe, easy to use, and the design is proven. The Secure II bed is used in Hospitals, Nursing Homes, Homes, and Teaching Schools. Secure II Beds can be configured to meet Hospital and Nursing Home needs. This unique bed is capable of being positioned for emergency situations and patient comfort. Teaching Schools get an easy to operate, safe, up to date hospital bed to train new nurses. The Secure II bed has optional features such as scale, bed exit, side rail communications, and IV Pole.

Safety is First with Stryker Secure II

Faced with increasing patient weight, age and acuity, clinicians unwittingly take on more risk as they are counted upon to provide superior care despite limited time and resources. Such expectations demand a patient handling solution that not only protects your patients, but also reduces the possibility of compromising your personal safety.

Easy to Use Proven Design: Stryker 3002 Secure II Bed Perfect for Varying Needs

Intuitive operation is paramount in today's fast-paced environment where time is critical. The Stryker 3002 Secure II bed features a wide range of clinical capabilities coupled with straight-forward, easy to understand controls make the Secure II an exceptionally convenient bed to use.

Fully customizable through a host of features designed around safety and ease of use, Secure II provides infinitely-adaptable functionality to meet your varying needs. Having a Stryker 3002 Secure II bed system that works for you allows you a greater opportunity to focus on your patient's recovery. Check out the other reconditioned Stryker products available from Piedmont Medical, Inc. including stretchers and new parts. Give us a call today! 1.800.433.3255 or 1.336.386.1040 Our helpful sales associates are waiting to hear from you.

- Overall Length: **93" (236.2 cm)**
- Overall Width (siderails up): **42 ½" (107.9 cm)**
- Overall Width (siderails down): **40" (102 cm)**
- Height Overall / High Position: **30" (76.2 cm)**
- Height Overall / Low Position: **16" (30.6 cm)**
- Caster Size: **6" (15.2 cm)**
- Fowler: **0° to 60°**
- Knee Gatch: **0° to 40°**
- Trendelenburg/Reverse Trendelenburg: **12°**
- Retraction: **14" (35.6 cm)**
- Patient Surface: **35" x 84" (88.9 cm x 213.4 cm)**
- Patient Weight: **500 lb. (226.8 kg)**

